

Multiplexed data analysis using QuPath

Nina Tubau Ribera

Centre for Dynamic Imaging – Advanced Technology and Biology

Walter and Eliza Hall Institute of Medical Research



Today's planning

- What is QuPath?
- Preprocessing: unmixing, removing autofluorescence

Live demo

- QuPath overview
- Pixel classifier
- Cell detection and classification
- Measurements
- Spatial analysis
- Scripting
- Advantages and drawbacks




What is QuPath?


- Open source for bioimage analysis
- Powerful annotation and visualisation tools
- Built-in algorithms for common tasks
- Interactive machine learning
- Compatibility
- Support for many file formats
- Scripting for deeper data interrogation





QuPath approach



Start with pixels
Billions of numbers




Identify objects
Cells, structures, regions...



Query the objects
Shape, size, number, intensity,
location, classification

Pre-processing

- Phenochart and InForm: Unmix fluorophores and autofluorescence



Let's get started !

Overview

Annotations Display Measurements

S | 0.67x | C | x,y

Project Image Annotations Hierarchy Workflow

Create project Open project Add images

image list

ProjectShirley/project.qpproj (2)

ER-4 TC_Scan1.qptiff - resolution #1

TNBC-3 TC_Scan1.qptiff - resolution #1

Image properties

Coordinates, pixel values
5485.80, 2251.49 µm
1, 0, 0, 2, 4, 0, 0, 25


800µm

Search entry in project

Pixel classifier

- Uses machine learning
- 1) Annotate a few areas
 - 2) Run pixel classifier


Classify > Pixel classification > Train pixel classifier



Cell detection


- Watershed

Analyze > Cell detection > Cell detection




Cell detection



- StarDist (run stardist.groovy)



Stardist(QuP
ath)



Watershed
(QuPath)




HALO

Cell classification

- Single measurement classifier



Classify > Object classification > Create single measurement classifier



Cell classification

- Composite classifier

Classify > Object classification > Create composite classifier




Define cell phenotypes




WEHI
brighter together

- Run `gui_phenotype.groovy`

Before running the script



After running the script



Types of measurements

GENERAL

Automatically measured

- Area, perimeter, centroid, intensities etc.

SPATIAL


Analyze>Spatial analysis> Select the measurement that you want to perform

- Distance to annotation 2D: distances between detection centroids and the closest annotation for each classification (eg. distance from every cell to closest area of tumour)
- Detect centroid distance 2D: distances between detection centroids for each classification. (eg. this may be used to identify the closest cell of a specified type)
- Delaunay cluster features 2D: Delaunay triangulation to detect objects based on their centroid locations. This helps identify clusters of objects neighbouring one another.
- Other possibilities (scripting needed): infiltration, nearest neighbours etc

Measurements


- Visualise measurements

Measure>Show measurement maps




Examples

- Distance to tumour edge



- Distance to specific type of cell



Batch processing and scripting

Possibility to run custom analysis or do batch processing

Workflow > Create Script

User forum, a good source of scripts: <https://forum.image.sc/tag/qupath>



Advantages and drawbacks

	QuPath	HALO
Advantages	Free and accessible Possibility to add any custom analysis	User-friendly
Drawbacks	Use of scripts	Limited to modules Only one license

Thank you!